SIEMENS

SONOCUR

	SP
Logbook	
Disassembly and Disposal tions	Instruc-
	© Siemens AG 1997 The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.


Register 3


Print No.: SPT3-130.163.01.01.02

Replaces: n.a.

English

Doc. Gen. Date: 04.97





Siemens AG

Fig. 1 Fig. 2




Fig. 3

General

The disposal of hazardous materials (materials that may be harmful to the environment) in a responsible manner is as important when disposing of the system/equipment as it is when replacing components and assemblies while the system is in use.

Disposal list

The following list shows the components/assemblies of the SONOCUR system and the hazardous materials which they contain.

The corresponding assemblies/parts are to be disposed of by an authorized waste disposal company.

Component/ Assembly	Location	Hazard/Hazardous material	Qty
ESWT cooling	1/Fig. 1	Coolant R 134 A	0.1 kg
system			
Board D3	4/Fig. 2	Lithium	1
Capacitor	2/Fig. 1	Oil	0.2 L
Spark gap	3/Fig. 1	Implosion	1
Spring	5/Fig. 3	Potential for injury	2

Tab. 1

TD SP 2 / Rottner SMS Iselin / Boyles

Disassembly and Disposal Instructions

This page intentionally left blank.

4

SONOCUR Register 3 SPT3-130.163.01 Page 4 of 4 Siemens AG Logbook Rev. 01 04.97 TD SP 1 Medical Engineering